


*The Gold Standard  
in EP Vessel Access*



**CLASSIC**SHEATH™  
Splittable Hemostatic Introducer System

# CLASSICSHEATH™

SPLITTABLE HEMOSTATIC INTRODUCER SYSTEM



## Product Features

**Robust valve design has been used for over 17 years**

**Has been shown to reduce back bleeding and risk of air embolism<sup>1</sup>**

**Color coded by French size for easy identification**

**Smooth dilator transition designed to reduce vein entry force**

**Snap-fit dilator secures assembly for introduction**

**Monolithic and atraumatic guide wire**

**Scored sheath for consistent peel away performance**

<sup>1</sup> Lee, H. Usefulness of a peelaway sheath introducer with a splittable hemostatic valve for placement of permanent pacemaker lead. *Catheterization and Cardiovascular Diagnosis*. 1997;40:326-327.

## CLASSICSHEATH Standard

Part Number	Sheath ID	French Size	Cap Color	Length
CLS-1006	0.083 in	6F	Green	13 cm
CLS-1007	0.092 in	7F	Orange	13 cm
CLS-1008	0.105 in	8F	Blue	13 cm
CLS-1008.5	0.112 in	8.5F	Light Blue	13 cm
CLS-1009	0.118 in	9F	Black	13 cm
CLS-1009.5	0.125 in	9.5F	White	13 cm
CLS-1010	0.131 in	10F	White	13 cm
CLS-1010.5	0.138 in	10.5F	White	13 cm
CLS-1011	0.144 in	11F	White	13 cm
CLS-1012	0.158 in	12F	White	13 cm
CLS-1012.5	0.164 in	12.5F	White	13 cm

## CLASSICSHEATH

Kit Contains:  
 ClassicSheath  
 Dilator  
 J-tipped Guide Wire  
 Syringe  
 18 Gauge Needle



## CLASSICSHEATH Long

Part Number	Sheath ID	French Size	Cap Color	Length
CLS-2507	0.092 in	7F	Orange	25 cm
CLS-2508	0.105 in	8F	Blue	25 cm
CLS-2509	0.118 in	9F	Black	25 cm
CLS-2510	0.131 in	10F	White	25 cm
CLS-2510.5	0.138 in	10.5F	White	25 cm
CLS-2511	0.144 in	11F	White	25 cm



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, MIDDLE EAST AND AFRICA (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 22

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone for specific country: Austria 0800 295 374 • Belgium 0800 72 906 (Dutch) 0800 73 172 (French) • Denmark 80 88 00 24 • France 0800 91 60 30

Finland 0800 770 586 • Germany 0800 182 0871 • Ireland (Republic) 1800 553 163 • Italy 800 897 005 • Luxembourg 8002 25 22 • Netherlands 0800 022 81 84

Norway 800 11629 • Sweden 020 792 445 • UK 0800 973 115

402876001/B ID 031913