




CARDIAC INTERVENTION

Hydrophilic Prelude SNAP™


Splittable Sheath Introducer

Now with
Hydrophilic
Coating

Requires
50%
less insertion
force



Easier to Break
requiring 66% less
break force than
its predecessor



Seamless
sheath-to-dilator
transition with
hydrophilic coating



Hydrophilic Prelude SNAP™

Splittable Sheath Introducer

The Hydrophilic Prelude SNAP features a **HYDROPHILIC** sheath coating that reduces insertion force by 50%.* A smooth sheath to dilator transition can help reduce vessel trauma and improve smooth navigation through torturous anatomy. The Hydrophilic Prelude SNAP contains the same **HIGH QUALITY** hub which requires 66% less break force than its predecessor making it **EFFORTLESS** to remove and helping ensure leads stay in place.

PRODUCT FEATURES

Robust Valve Design

Tight seal to protect against air embolism and blood loss.

Improved Grip

Expanded wings with non-slip, ridged surfaces facilitate a tighter, more assured grip.




Easier to Break

Requires 66% less break force than its predecessor.*

Seamless Transition

Smoother sheath-to-dilator transition to reduce vessel trauma and promote ease of insertion.



Competitors' sheath-to-dilator transitions

Our seamless sheath-to-dilator transition

Hydrophilic Coating

Covers 75% of distal sheath length for ease of insertion while maintaining stability.

KIT CONTENTS

- Splittable Sheath
- Dilator
- Guide Wire 0.038" x 50 cm or 0.038" x 80 cm
- 18G Needle
- 12 mL Syringe

ORDERING INFORMATION

Catalog Number with Sideport	French Size	Sheath Length	Hub Color
PLSH-1006	6	13 cm	Green
PLSH-1007	7	13 cm	Orange
PLSH-1008	8	13 cm	Blue
PLSH-1008.5	8.5	13 cm	Light Blue
PLSH-1009	9	13 cm	Black
PLSH-1009.5	9.5	13 cm	Gray
PLSH-1010	10	13 cm	Fuchsia
PLSH-1010.5	10.5	13 cm	Light Fuchsia
PLSH-1011	11	13 cm	Yellow
PLSH-1012	12	13 cm	Brown
PLSH-1012.5	12.5	13 cm	Light Brown
PLSH-2506	6	25 cm	Green
PLSH-2507	7	25 cm	Orange
PLSH-2508	8	25 cm	Blue
PLSH-2509	9	25 cm	Black
PLSH-2510	10	25 cm	Fuchsia
PLSH-2510.5	10.5	25 cm	Light Fuchsia
PLSH-2511	11	25 cm	Yellow

Catalog Number without Sideport	French Size	Sheath Length	Hub Color
PLSHX-1006	6	13 cm	Green
PLSHX-1007	7	13 cm	Orange
PLSHX-1008	8	13 cm	Blue
PLSHX-1008.5	8.5	13 cm	Light Blue
PLSHX-1009	9	13 cm	Black
PLSHX-1009.5	9.5	13 cm	Gray
PLSHX-1010	10	13 cm	Fuchsia
PLSHX-1010.5	10.5	13 cm	Light Fuchsia
PLSHX-1011	11	13 cm	Yellow
PLSHX-1012	12	13 cm	Brown
PLSHX-1012.5	12.5	13 cm	Light Brown
PLSHX-2506	6	25 cm	Green
PLSHX-2507	7	25 cm	Orange
PLSHX-2508	8	25 cm	Blue
PLSHX-2509	9	25 cm	Black
PLSHX-2510	10	25 cm	Fuchsia
PLSHX-2510.5	10.5	25 cm	Light Fuchsia
PLSHX-2511	11	25 cm	Yellow

*Data on file.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Merit.com

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Sweden
020 792 445

Switzerland
(Deutsch)
+41 225180252
(Français)
+41 225948000
(Italiano)
+41 225180035
UK
0800 973 115