



### Functions:

- The tension of the materials can be used to open the incision & better operation vision can be obtained
- The Soft tissue retractor is closely combined with incision
- Latex-free silicone rubber & polyurethane offers ultimate versatility goods


### Benefits:

- Enhances working space, improves visibility and exposure
- Distributes force evenly, potentially minimizing tissue trauma & pain
- Reduces the need for hand-held retractors, maximizing enhances surgical efficiency
- Allows visualization and protects wound margins tissue
- Maintains moisture at the incision site and reduces tissue damages to the minimum
- Significantly reduces the risk of incision infection
- Effectively prevent the contact between cut organs and abdomen cavity
- Reduces risk of embolic fat particles in the instruments, sutures and knot pushers
- Very simple, convenience, practicality facilitates faster and effortless placement
- Height flexibility and changeably, accommodates a wide range of patients
- Disposable, ready to use, safe



## Soft Tissue Retractor With Fixed Height

### Fixed Height




Type	A(mm)	B(mm)	C(mm)	D(mm)
HK-50/40-15/20 ( D )	50	40	15	20
HK-50/40-25/25 ( D )	50	40	25	25
HK-60/50-25/25 ( D )	60	50	25	25
HK-70/70-35/25 ( D )	70	70	35	25
HK-100/80-60/25 ( D )	100	80	60	25
HK-110/100-60/25 ( D )	110	100	60	25
HK-120/100-60/25 ( D )	120	100	60	25




HTKD Medical



*Atraumatic retraction  
Tightness attach*



*Simple operate  
Superior exposure*



**Soft Tissue Retractor With Changeable Height**



**Changeable Height**

Type	A <sup>∅</sup> (mm)	B <sup>∅</sup> (mm)	C <sup>∅</sup> (mm)	D <sup>∅</sup> (mm)
HK-280/290-280/100 ( B )	280	290	280	40-150
HK-300/310-300/100 ( B )	300	310	300	40-150
HK-320/330-320/100 ( B )	320	330	320	40-150

