

# BUILT FOR STRENGTH

**-HERO<sup>®</sup>**  
**GRAFT**  
**Venous Outflow Component**

**HERO1001VOC**

- Integrated, crush resistant radiopaque marker band
- Tougher blend of silicone for improved pushability



## HeRO Graft

HeRO Graft (Hemodialysis Reliable Outflow) is a fully subcutaneous AV access solution clinically proven to maintain long-term access for haemodialysis patients with **central venous stenosis**.

## Silicone Venous Outflow Component

- No venous anastomosis
- Reinforced: 48 braid nitinol for kink & crush resistance
- Removable and replaceable
- Radiopaque band (at distal tip)

## Treatment Algorithm

Failing AVF or AVG due to central venous stenosis



Catheter-dependent patients

Compared to the HERO1001

**the HERO1001VOC has a:**

- Tougher blend of silicone for improved pushability
- Nitinol radiopaque marker band which is:
  - Integrated
  - Crush resistant



## Ordering Information

Product Code	Component	Diameter (ID)	Length
<b>HERO1001VOC</b>	Venous Outflow Component (with nitinol marker band)	5mm	40cm (customisable)
<b>HERO1002</b>	Arterial Graft Component	6mm (ePTFE); 6mm - 5mm (connector)	53cm (connector: 3cm)
<b>HERO1003</b>	Accessory Component Kit	N/A	N/A

**Before using: Refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.**



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

[merit.com](http://merit.com)

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Singapore  
+65 67536311

South Korea  
+8210 2205 0707

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115